DAYLIGHT GUARANTEE/WARRANTY

This product is guaranteed against mechanical and electrical defects. For full Daylight guarantee Terms & Conditions please refer to the Daylight website. Faulty products should be returned to the place of purchase and proof of purchase must be presented. Please quote the model number in all correspondence. Returned products damaged by poor packaging are the responsibility of the sender. During the guarantee period, Daylight will repair or replace the product free of charge provided the fault was not caused by misuse, improper handling, negligence, normal wear and tear, accidental damage, or by any modifications made by the purchaser. The guarantee does not cover consumables such as bulbs, tubes or batteries. This guarantee does not affect your statutory rights.

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


SUOMALAINEN

DANSK


EN - This light source (LED only) is replaceable by manufacturer, by service personnel authorised by the manufacturer or by a comparably qualified person.

This product contains a light source (LED only) of energy efficiency class < D >

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EN - The Control gear is replaceable by the end user.

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A GUIDE TO LIGHT THERAPY

Light therapy can boost your mood, deepen your sleep, energise you and sharpen your concentration.

You might start to notice that shorter days and reduced sunlight are affecting your energy and mood. This is often referred to as winter blues, starting in autumn and lasting several weeks. For some individuals it can last throughout the winter.

If you're feeling blue or experiencing SAD, regular light therapy can alleviate your symptoms, helping you to feel back to your usual self.

During spring and summer, as natural light levels increase, you may wish to reduce the regularity of your sessions, though if you feel your symptoms coming back, you should adjust your light therapy schedule accordingly.

The key factors to consider when using light therapy are Light Intensity, Distance & Time

To simulate the brightness of natural sunlight, we recommend using the SAD light on full brightness at arm's length (50-70cm), ideally every day for minimum 30-60 minutes.

Example session times and Lux output:

30 minutes with the lamp approximately 16cm (6") away from you (10,000 Lux)

60 minutes with the lamp approximately 24cm (9") away from you (5,000 Lux)


120 minutes with the lamp approximately 36cm (14") away from you (2,500 Lux)

With the light closer to you, a greater intensity of light reaches your eyes. However, as individual sensitivity to light varies, we recommend finding the right distance and brightness that works best for you. If you find the light too bright, simply place the lamp further away from you while increasing your session time.

Alternatively, use the dimmer to reduce the brightness and increase the session duration to ensure you're receiving the right amount of light.

Positioning of the light is important to get the best results and to ensure optimal comfort

Position the light to your side, with the shade directed towards you. Be sure to keep your eyes open as the treatment will not work on closed eyes. Do not stare into the light, it is sufficient that the light reaches your eyes indirectly.

Light therapy is at its most effective when used first thing in the morning.

Avoid using your lamp in light therapy mode in the evening as you may have difficulty falling asleep. Your doctor can help you find the most beneficial light therapy schedule for you.

This is not a medical guide - please see your doctor in relation to any medical conditions and follow their guidelines. If you do not see an improvement in your symptoms after a few weeks, please seek advice from your doctor.

SET-UP AND OPERATION

1. Remove all packaging from around the lamp
2. The operator should insert the power adapter jack into the socket at the rear of the lamp and plug the adapter into a conveniently located electrical socket.
3. Press button (1) to switch ON and press again to switch OFF.
4. Rotating button (1) clockwise will increase brightness, and anti-clockwise will decrease brightness.
5. The colour of the light can be changed by touching the colour change icon (2).
For every touch, the light colour will cycle between daylight (6,000K), cool light (4,000K) and warm light (2,700K).
6. To start 'Light Therapy Mode' touch the sun icon (3).
The lamp will flash once to indicate the start of light therapy and will automatically switch to daylight (6,000K) and on maximum brightness for 30 minutes.
7. After 30 minutes of 'Light Therapy Mode' the lamp automatically turns off. If you would like to continue your light therapy session, simply touch the sun icon (3) and the lamp starts another 30 minute session. Alternatively,

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pressing the button (1) switches the lamp back into 'Desk Mode', brightness and light colour can be adjusted as described above.

8. In both Desk or Therapy Modes the shade and arm can be adjusted, angled and rotated as required to enable optimum position of the light. When in Desk Mode the timer function will not operate and the lamp will remain on until turned off.

9. The base of the lamp has a USB port (4) with 5V / 1.2A output. It is suitable for charging mobile phones and other similar USB powered devices.

Simply insert the USB cable that is supplied with your device into the USB port. The port provides power whether the lamp is switched on or off.

CAUTION: DO NOT CONNECT USB DEVICES WHICH REQUIRE MORE POWER THAN 5V / 1.2A

Tip 1: As with 'Desk Mode', while in 'Light Therapy Mode' the brightness can be reduced by rotating the button (1).

Tip 2: During 'Therapy Mode' if warmer light colour is preferred, the colour change icon (2) can be used to select cool light or warm light. Adjusting the brightness or light colour does not affect the 30 minute timer, the timer is started from the moment the sun icon (3) is touched. Retouching the sun icon at anytime within the 30 minutes will restart the 30 minute timer and default back to daylight maximum brightness.

Please note: Both cool light and warm light provide a good therapy light and some individuals prefer the softer warmer light, however, for optimum results in the shortest time we recommend daylight is used for light therapy sessions.

PRECAUTIONS

Light Therapy is a safe and effective treatment for winter blues such as Seasonal Affective Disorder (SAD) as well as for helping relieve the symptoms of poor sleep / wake patterns and low energy levels.

Before using the lamp we advise seeking advice from your doctor if any of the following apply to you:

- have been advised to avoid bright light because of a particular medical condition
- taking medication that may cause photosensitivity
- had recent eye surgery or have had previous eye problems
- have been diagnosed with a mood disorder, depression or sleep problems
- take particular care if you have been prescribed any medication for any of the above, as light therapy can also affect your mood and affect how your medication works

Light Therapy can have some side effects, but these are usually mild and can be resolved by using the light further away / lowering the brightness.

For some individuals a brief burst of bright light for 30min is very beneficial and is preferred. For others a lower intensity light, used for longer is more comfortable and gives better results. Side effects that might be experienced include: eye strain, headache and nausea. To minimise the possibility of side effects do not stare directly into the light, start by gradually building up the time you spend in front of the light and consider a lower light level and greater distance to give your eyes and body time to adapt. It is important to listen to your own body and find the brightness, distance and treatment time that works best for you.

This is not a medical guide or directions – If your doctor has recommended light therapy always seek advice from them and follow their instructions. If you have any concerns about anything raised in this guide please speak to your doctor.

CLEANING

BEFORE CLEANING MAKE SURE THE LAMP IS DISCONNECTED FROM THE ELECTRICAL SUPPLY.

If the lamp requires cleaning simply wipe with a damp cloth. Do not use large amounts of liquid or spray cleaners as they may enter the lamp and affect electrical safety.

MAINTENANCE

No maintenance is required. The LEDs are not user or service replaceable as they designed to last for the lifetime of the product.

ELECTRICAL SAFETY

FOR INDOOR USE ONLY

CAUTION: If the cable of this luminaire becomes damaged, it should be exclusively replaced by a suitably qualified person in order to avoid hazard.

IN CASE OF DOUBT CONSULT A QUALIFIED ELECTRICIAN

FRANÇAIS

IMPORTANT SAFETY INSTRUCTIONS – PLEASE SAVE FOR FUTURE REFERENCE

THANK YOU FOR PURCHASING THIS DAYLIGHT PRODUCT. TO ENSURE CORRECT AND SAFE OPERATION PLEASE READ THE FOLLOWING AND RETAIN FOR FUTURE REFERENCE.

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